

1. Identification**Product name****HigherPurity™ Blood DNA Extraction Kit**

(Isolation of DNA from 200 µl-2ml Whole Blood)

Cat. No**AN0043-M (50 - 5 reactions)****Cat. No****AN0043 (150 - 15 reactions)****Cat. No****AN0043-XL (250 - 25 reactions)****2. Description**

HigherPurity™ Blood DNA Extraction kit is a simple and rapid method for high-quality genomic DNA purification from various sources, including: whole blood, buffy coat* and cultured cells*.

The procedure includes: lysis, protein removal, DNA precipitation, washing and hydration.

*Extraction protocols are available at the Canvax Biotech website.

3. Kit Components

Item	Quantity		
	AN0043-M	AN0043	AN0043-XL
S1 Buffer	40 mL	2 x 60 mL	3 x 60 mL
S2 Buffer	12 mL	36 mL	60 mL
S3 Buffer	5 mL	12 mL	20 mL
Proteinase K*	11 mg	30 mg	2 x 30 mg
EB Buffer	5 mL	12 mL	20 mL

*Dissolve Proteinase K in water to obtain a 20 mg/mL stock solution. **The Proteinase K solution can be stored for several days at 2–8 °C.** For longer-term storage, the unused portion of the solution may be stored in aliquots at –20 °C until needed.

4. Storage

Store the kit at room temperature. If any kit reagent forms a precipitate, warm at 55–65 °C until the precipitate dissolves and allow to cool to room temperature before use.

5. Features

- **Safe:** no phenol-chloroform extraction.
- **Efficient:** 3–6 µg of genomic DNA from a 200 µl blood sample (30–90 µg of genomic DNA from a 2ml).
- **Ready to use** genomic DNA, in all molecular biology applications.

6. Applications

High molecular weight genomic DNA purified with the kit is suitable for direct use in all common molecular biology applications: PCR, cloning, DNA sequencing, Southern blot analysis, etc.

7. Quality Certifications

HigherPurity™ Blood DNA Extraction kit is tested on a lot-to-lot basis by isolating total DNA from 200µl of whole human blood. DNA purified is analyzed by:

- Spectrophotometer: Ratio 260/ 280 (1.6–1.8)
- Agarose gel electrophoresis.



8. Further information

Product Use Limitations	This product is developed, designed, and sold exclusively only for research purposes use. The product was not tested for use in diagnostics or for drug development, nor is it suitable for administration to humans or animals.
Safety Information	When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles. For more information, please consult the appropriate material safety data sheets (MSDS). These are available online in convenient and compact PDF format at www.canvaxbio.com where you can find, view, and print the MSDS for each CANVAX kit.

DETAILED PROTOCOL

Choose ■ if processing 200 µl blood samples // Choose ▲ if processing 2 ml blood samples.

- Transfer up to **■ 200 µL / ▲ 2ml** of sample to a **■ 1.5 ml / ▲ 15ml tube**.
- Add **■ 600 µL / ▲ 6ml buffer S1** to the sample and vortex gently or invert tube 6-8 times and leave to incubate for **2-3 minutes / 5-8 minutes** at room temperature
- Centrifuge at **■ 13000 rpm / ▲ 2500 rpm** for **1-2 minutes**. Remove the supernatant using pipette and avoiding damaging the cell visible pellet and leaving **■ 50 µl / ▲ 200 µl** of residual liquid.
- Vortex the tube vigorously until the white blood cells are resuspended (10-15 seconds). **This process will help to optimize the cell lysis in the following step.**
- Add **■ 5-10 µL / ▲ 50-100 µl** proteinase K and **■ 200 µL / ▲ 2ml buffer S2** and mixing by pipetting. Transfer all the mix into a new **■ 1.5 ml / ▲ 15ml tube**.
- Incubate in a water bath at 55 °C for 0.5-1 hour, and then cool to room temperature.
- Add **■ 67 µL / ▲ 670 µl buffer S3** and mixing with vortex vigorously for 20 seconds.
- Centrifuge at **■ 13000 rpm for 1 minute / ▲ 2500 rpm for 5 minutes**. A dark brown pellet should be visible. If no pellet is observed, incubate on ice for 5 minutes and centrifuge again.
- Transfer the supernatant to a new **■ 1.5 ml / ▲ 15ml tube** containing **■ 200 µL / ▲ 2ml** isopropanol. Mix by gentle inversion 50 times.
- Centrifuge at **■ 13000 rpm for 1 minute / ▲ 2500 rpm for 3 minutes** and remove the supernatant. The DNA will be visible as a small white pellet.
- Wash with **■ 200 µL / ▲ 2ml 70% ethanol** and centrifuge at **■ 13000 rpm for 1 minute / ▲ 2500 rpm for 1 minutes**.
- Remove the supernatant using a pipette and dry the pellet with the tube inverted on absorbent paper for 5-10 minutes.
- Resuspend the DNA adding **■ 65 µL / ▲ 200µl of Buffer EB**.
- Store DNA at -20 ° C.

