

1. Identification

Product name	WideUSE™ Plasmid Purification MiniSpin Kit
Cat. No	AN0068 (50 reactions)
Cat. No	AN0069 (100 reactions)
Cat. No	AN0069-XL (250 reactions)

2. Description

The WideUSE™ Plasmid Purification Kit is designed for rapid and cost-effective small-scale preparation of high-quality plasmid DNA from recombinant E. Coli cultures. By combining silica-binding technology and the convenience of a spin column format, up to 24 µg of high copy plasmid DNA can be recovered from 1-5 mL of E. coli culture in less than 30 minutes. The extraction process includes an initial SDS/alkaline lysis step with the appropriate Buffer to liberate plasmid DNA from the cell followed by adsorption of the DNA onto silica in the presence of high salts. Contaminants are then removed by a spin-wash step. Finally, the bound DNA is eluted in water or Tris-EDTA buffer.

3. Kit Components

Item	Quantity		
	AN0068 (50 rxn)	AN0069 (100 rxn)	AN0069-XL (250 rxn)
WideUSE minispin columns	50	100	250
Collection tubes (2 mL)	50	100	250
S-I Buffer	8 ml	16 ml	40 ml
S-II Buffer	8 ml	16 ml	40 ml
S-III Buffer	8 ml	16 ml	40 ml
Binding Buffer	30 ml	60 ml	150 ml
Washing Buffer*	8 ml	16 ml	40 ml
Elution buffer	6 ml	10 ml	20 ml
RNase A	80 µl	160 µl	400 µl

*Add the volume ethanol (96%-100%) specified [Not included] to WB Buffer prior to initial use (see the label on the bottle for a volume indication). After ethanol has been added, mark the bottle to indicate that this step has been completed.

4. Features

- **High yields** of up to 24 µg of DNA suitable for all molecular biology procedures.
- No phenol-chloroform extraction.
- **Ready to use** plasmid DNA.
- **Just a few minutes** procedure.
- **Mini format**



DATA SHEET

Version: 03
Revision date: 11/04/2023

Canvax Reagents, S.L.U.
Luis de Mercado Street, 19
Boecillo Technological Park
47151, Boecillo
Valladolid, Spain.

Tlf: +34 983 54 85 63
info@canvaxbiotech.com

www.canvaxbiotech.com

5. Storage specifications

WideUSE™ Plasmid Purification Kit should be stored at room temperature (15–25°C) for up to 24 months without any reduction in performance. Store the RNase A at -20°C. After addition of RNase A to Buffer S-I reagent can be stored at 4°C.

6. Applications

All molecular biology applications, such as:

- Digestion with restriction enzymes.
- Automated sequencing.
- PCR template.
- Bacterial transformation.

7. Quality Certifications

Plasmid Purification Kit is tested for the isolation of any plasmid DNA from transformed E.coli. The quality of purified DNA is analysed by:

- Ratio 260/ 280.
- Agarose gel electrophoresis.
- Digestion with restriction endonucleases.

8. PRODUCT USE LIMITATION

Product Use Limitations This product is developed, designed, and sold exclusively only for research purposes use. The product was not tested for use in diagnostics or for drug development, nor is it suitable for administration to humans or animals.

Safety Information When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles. For more information, please consult the appropriate material safety data sheets (MSDS). These are available online in convenient and compact PDF format at www.canvaxbio.com where you can find, view, and print the MSDS for each CANVAX kit.



DATA SHEET

Version: 03
Revision date: 11/04/2023

Canvax Reagents, S.L.U.
Luis de Mercado Street, 19
Boecillo Technological Park
47151, Boecillo
Valladolid, Spain.

Tlf: +34 983 54 85 63
info@canvaxbiotech.com

www.canvaxbiotech.com

PROTOCOL

PREPARATION OF WORKING SOLUTIONS

Before starting the protocol prepare the following reagent:

- ✓ To prepare the S-I Buffer add the vial of **RNAse A**.
- ✓ Add the volume ethanol (96%-100%) specified [Not included] to Washing Buffer prior to initial use. After ethanol has been added mark the bottle to indicate that this step has been completed.

ASSAY PROCEDURE

1. Pour the culture in a 1.5 ml centrifuge tube and harvest the bacterial cells by centrifugation
at 13000 rpm for 2 minutes. For low-copy plasmids should be collect 3 mL of culture and
using 2 volumes of each solution to obtain good yields.
2. Resuspend the bacterial pellet in 100 μ L of **Buffer S-I**
3. Add 100 μ L of **Buffer S-II**, mix thoroughly by inverting the tube 6 times.
4. Add 100 μ L of **Buffer S-III**, mix thoroughly by inverting the tube 8 times.
5. Centrifuge at 13000 rpm in a microcentrifuge for 10 min. Recover supernatant containing
plasmid DNA promptly into a 1.5 ml centrifuge tube.
6. Add 500 μ L of **Binding Buffer**, mix by inverting the tube several times.
Incubate at room
temperature for 5 min.
7. Apply the supernatants from step 6 to the **WideUSE spin column** by
decanting or
pipetting.
8. Centrifuge at 5500 rpm (9500g) for 90 seconds. Discard the flow-through.
9. Wash the **WideUSE spin column** by adding 700 μ L **washing Buffer** and
centrifuging at
5500 rpm for 90 s. Discard the flow-through.
10. Place the **WideUSE spin column** in a collection tube and add 500 μ L of
Isopropanol pure.
11. Centrifuge at 5500 rpm for 90 s. Discard the flow-through.
12. Again Centrifuge at 13000 rpm for 90 s. This step helps to remove traces of
isopropanol.
13. Place the **WideUSE spin column** into a new, labelled 1.5 microcentrifuge
tube and pipet
50-60 μ L **Elution buffer** directly into the membrane or pre-warm water. Close
the cap and
incubate for 1 minute at room temperature.
14. Centrifuge at 13000 rpm for 1 minute to elute DNA.

