

DATA SHEET

Version: 2
Revision date: 17/02/2023

Canvax Reagents, S.L.U.
Luis de Mercado Street, 19
Boecillo Technological Park
47151, Boecillo
Valladolid, Spain.

Tlf: +34 983 54 85 63
info@canvaxbiotech.com

www.canvaxbiotech.com

1. Identification

Product name **Tris Buffered Saline with Tween™ 20 (TBS-T), pH7.6**
8L (16 tablets)

Cat. No **BR0044**

2. Description

TBS Tween™20 (TBS-T) is the most commonly used wash buffers or blocking buffer used for applications such as Western blot and ELISA assays. TBST is an optimal formulation of pH stabilizers, salts and detergents designed to remove excess material, decreasing non-specific background staining. It is preferably used with alkaline phosphatase or peroxidase-conjugated antibodies.

3. Specifications

Chemicals: Analytical grade.
Format: Exactly pre-weighed tablets
Composition: : 0.15 M NaCl, 0.050 M Tris-HCl, 0.05% Tween™ 20
pH: 7.6 ± 0.2 at 25 °C

4. Shipping and Storage specifications

TBS-T buffer is shipped at room temperature. Store the product in a dry place at room temperature.

5. Applications

- Washing nitrocellulose membrane in Western blot and microtiter plate wells in ELISA assays
- Blocking buffer for plate-based assays

6. Directions for use

Dissolve one tablet in 500 ml of water.

7. Further information

Product Use Limitation This product is developed, designed and sold exclusively only for research purposes use. The product was not tested for use in diagnostics or for drug development, nor is it suitable for administration to humans or animals.

Disclaimer The information provided in this Data Sheet is correct to the best of our knowledge and belief at the date of publication. This information is intended only as a guide and should not be taken as a warranty or quality specification. Canvax Reagents S.L.U. shall not be held liable for any damage resulting from handling or from contact with the above product.

