

DATA SHEET

Version: 03


Revision date: 15/06/2023

1. Identification

Product name	pSpark® III (20ng/μL) 20 reactions
Cat No	C0003

2. Description

pSpark® III is a highly efficient, accurate and easy-to-use DNA cloning system that combines Ampicillin and Kanamycin resistance. Ideal for cloning PCR products amplified from any plasmid vector without the need to gel-purify bands to eliminate the background due to the template vector used for PCR.



3. Composition

Item	Quantity
pSpark® III DNA cloning vector	20 rxn (20ng/ μL)
10X PEG 6000 solution	150 μL
1kb control insert	5 μL
T4 DNA ligase (5 Weiss U/μL)	20 μL
10x T4 DNA ligase Buffer	100 μL

4. Features

- **Unprecedented high cloning efficiency:** > 2,500 positive colonies expected under optimal conditions.
- **Time-saving protocol:** no hidden steps such as phosphorylation, just ligation after PCR and transformation.
- **Powerful:** obtain 5x more positive colonies using 10x less DNA insert.
- **Easy-to-use:** eliminate recombinant screening due to its < 1% background, avoiding "suicide" strategies from toxic genes.
- **High stability:** eliminates cloning bias or pitfalls.
- **Great versatility:** compatible with any protocol, proofreading polymerase, competent cells, ligation time or primers.
- **Sensitive:** clone from 50 bp insert to up to 14 kb with just 5ng per kb of insert.
- **High cost-saving:** reduces your cloning costs as no expensive phosphorylated primers are needed.
- **Eliminates positive selection vector.**



DATA SHEET

Version: 03

Revision date: 15/06/2023

Tlf: +34 983 54 85 63
info@canvaxbiotech.comwww.canvaxbiotech.com

5. Storage specifications

Store in a non-frost-free freezer at **-20°C**.

6. Applications

- Cloning directly from PCR using plasmid cloned genes as template.
- Unpurified PCR cloning.
- Cloning of high-fidelity PCR amplified products.
- Production of ssDNA.
- Blue/white screening for recombinants.
- In vitro transcription from T7/SP6 dual-opposed promoters.

7. Recommended protocol

1. Amplification

Amplify DNA with any proofreading polymerase following suggested protocol. Gel purifies the insert or check by electrophoresis the quality of amplified DNA. If a single band appears on the gel and template DNA is not a vector with Ampicillin resistance, like the pSpark® III DNA cloning vector, then gel purification is not needed.

2. Ligation

Suggested ratio of insert to vector is 5:1. Use the equation **ng insert = 33ng x kb** (of insert, that is, for a 1kb insert only 33 ng are needed) or use **1 µL** of unpurified PCR product for directly ligation.

Mix:

1 µL of pSpark® III DNA cloning vector (20ng/µL)

1 µL 10x T4 DNA ligase buffer

1 µL 10x PEG 6000 solution

X µL of insert

1 µL of T4 DNA ligase (5 Weiss units)

H₂O to **10 µL**

Incubate the reaction 60 minutes at 22°C (a ligation time from 10min to overnight could be used)

3. Transformation

Transformation can be made with your regular protocol, competent cells and selective media. As pSpark® III DNA cloning vector has Ampicillin resistance, this antibiotic is needed. Also, as pSpark® III DNA cloning vector has blue/white screening capability, IPTG/XGal plates are recommended.

A suitable transformation protocol is:

- Thaw competent cells on ice, take **50 µL** and put into a cool 1.5 mL tube.
- Add **10 µL of the ligation mixture**, mix by gently flicking with fingertips and incubate on ice for 30 minutes.
- Heat-shock the cells at **exactly 42°C for exactly 45 seconds** (do not shake) and then transfer to ice for 2 minutes.
- Plate the cells into media with antibiotic, IPTG and XGal and incubate overnight at 37°C.



DATA SHEET

Version: 03

Revision date: 15/06/2023

Tlf: +34 983 54 85 63
info@canvaxbiotech.comwww.canvaxbiotech.com

4. Analysis of transformants

Either colony PCR or plasmid purification and digestion could be used for analysis of transformants. pSpark® III DNA cloning vector has pUC/M13 forward and reverse sequencing primers binding site for sequencing or amplification of insert. Insert size when amplified with pUC/M13 fw and rev primers is about 200bp longer than the real size of insert.

5. Expected results

Up to 2500 positive white colonies for a 60 min ligation of 33 ng of supplied control insert (1 kb) and less than 20 negative blue colonies when using cells with a transformation efficiency of 4×10^7 cfu/ug.

7. Further information

Product Use Limitations	This product is developed, designed and sold exclusively only for research purposes use. The product was not tested for use in diagnostics or for drug development, nor is it suitable for administration to humans or animals.
Disclaimer	The information provided in this Data Sheet is correct to the best of our knowledge and belief at the date of publication. This information is intended only as a guide and should not be taken as a warranty or quality specification. Canvax Reagents S.L.U. shall not be held liable for any damage resulting from handling or from contact with the above product.

