

# DATA SHEET

Version: 03  
Revision date: 27/03/2023

## 1. Identification

Product name	<b>GreenSafe DNA gel stain</b>
	500 µL
Cat. No	<b>E0206</b>

## 2. Description

**GreenSafe DNA gel stain** product replace toxic Ethidium Bromide (EtBr) for the visualization of doublestranded DNA, single-stranded DNA, and RNA in agarose and polyacrylamide gel electrophoresis. GreenSafe DNA gel stain has 2 fluorescence excitation wavelengths in the UV range (~270nm; ~290nm) and one in the blue light range (~485nm). Maximum fluorescence emission is at ~525nm (green). Therefore, GreenSafe DNA gel stain is compatible with a large variety of gel documentation systems.

## 3. Composition

Item	Quantity
GreenSafe DNA gel stain	500uL

## 4. Features

- **Easy to Use:** you can directly replace EB with without changing your existing imaging system.
- **Safe:** Non-carcinogenic by the AMES test.
- **Sensitive:** Increase your sensitivity by reducing nonspecific background fluorescence.

## 5. Storage specifications

You may store the **GreenSafe DNA gel stain** at any temperature between 2°C to 8°C protected from light. Do not freeze!

## 6. Applications

Non-carcinogenic alternative to Ethidium bromide.

## 7. Pre-casting protocol

**500 µl of this stain is sufficient for 10 L of agarose gel.**

1. Prepare 100 ml of agarose gel solution (concentration from 0.8 - 3.0 %) and heat until the solution completely clear and no small floatin particles are visible.
2. Add 4-6 µl of the stain to the gel solution and mix it gently.
3. Cool the gel to 50 – 60 °C and cast the gel into the gel tray.
4. When the gel is solid, load the samples and perform electrophoresis.
5. Detect the bands under UV illuminator. Can also be viewed with Blue LED light. Yellow or green gelatin or cellophane filters should be used for photography.



## DATA SHEET

Version: 03

Revision date: 27/03/2023

### 8. Poststaining

1. For <0.5 cm thick agarose gels, 10-15  $\mu$ L of stain should be used per 100 mL of buffer
2. Incubate the gel in staining solution for 10 – 30 minutes. Staining time varies with the thickness of the gel and percentage of the agarose.
3. The post-staining solution may be used 2-3 times. Staining solution that is going to be reused should be preferably stored at room temperature in the dark

### 9. Further information

<b>Product</b>	This product is developed, designed and sold exclusively only for research purposes use.
<b>Use</b>	The product was not tested for use in diagnostics or for drug development, nor is it suitable
<b>Limitations</b>	for administration to humans or animals.
<b>Disclaimer</b>	The information provided in this Data Sheet is correct to the best of our knowledge and belief at the date of publication. This information is intended only as a guide and should not be taken as a warranty or quality specification. Canvax Reagents S.L.U. shall not be held liable for any damage resulting from handling or from contact with the above product.

