

# DATA SHEET

Version: 02  
Revision date: 13/03/2023Tlf: +34 983 54 85 63  
[info@canvaxbiotech.com](mailto:info@canvaxbiotech.com)[www.canvaxbiotech.com](http://www.canvaxbiotech.com)

## 1. Identification

**Product name****T7 Digoxigenin RNA Labeling Kit**

RNA Synthesis with Modified Nucleotides

30 reactions

**Cat. No.**


RS024

## 2. Description

**Digoxigenin (DIG)** is a steroid hapten derived from plants of the genus *Digitalis* (*D. lanata*, and the purple foxglove *D. purpurea*). It is widely used for generating non-radioactive probes and is an alternative to biotinylation. Like biotin, DIG-labeled probes support both chromogenic and chemiluminescence detection formats.

RNA probes are synthesized by in vitro transcription with **DIG-11-UTP**. These nucleotides are linked via a spacer arm to DIG. The resulting DIG-labeled molecules then function as hybridization probes in much the same manner as any other type of probe.

**T7 Digoxigenin RNA Labeling Kit** is designed to generate randomly Digoxigenin-modified RNA probes by in vitro transcription with Digoxigenin-11-UTP



## 3. Kit components

Item	Volume
T7 RNA Polymerase mix *	2x 40 µl
T7 10X Reaction Buffer **	200 µl
ATP (100 mM)	100 µl
GTP (100 mM)	100 µl
CTP (100 mM)	100 µl
UTP (100 mM)	100 µl
DIG-11-UTP (10mM)	25 µl
DTT (100 mM)	100 µl
Control template 1 (200 ng/µl)	10 µl
Nuclease-Free Water	1.2 ml

\*Include RNase inhibitor

\*\*HEPES-based

## 4. Recommended Storage Condition

Store all kit components at -20 °C. Avoid freeze/thaw cycles.

**T7 Digoxigenin RNA Labeling Kit** components are stable up to the expiry date specified with the product.

## DATA SHEET

Version: 02

Revision date: 13/03/2023

### 5. Applications

- In situ hybridization
- Southern blots
- Microarray hybridization
- Chromosome mapping

### 6. Further information

- Product Use** This product is developed, designed and sold exclusively only for research purposes use. The product was not tested for use in diagnostics or for drug development, nor is it suitable for administration to humans or animals.
- Limitations**
- Disclaimer** The information provided in this Data Sheet is correct to the best of our knowledge and belief at the date of publication. This information is intended only as a guide and should not be taken as a warranty or quality specification. Canvax Reagents S.L.U. shall not be held liable for any damage resulting from handling or from contact with the above product.



## DATA SHEET

Version: 02

Revision date: 13/03/2023

### IN VITRO TRANSCRIPTION PROTOCOL

#### Materials to Be Supplied by the User:

- T7 Promotor-containing DNA template
- RNA purification tools
- RNase-free DNase I

#### Important Notes before starting:

- Take all precautions to avoid RNase contamination during the whole process. We strongly recommend wearing gloves and using nuclease-free tubes and reagents to avoid RNase contamination.
- Reactions are typically 20 µL but can be scaled up as needed.
- All reagents should be spin briefly before opening to prevent loss and/or contamination.
- The quality and purity of input RNA and resulting cDNA determine the yield of DOG-labeled RNA.
- The recommended molar ratio of T7 Digoxigenin RNA Labeling Kit to standard NTP is 1:3 or 1:2.
- The kit contains sufficient reagents for 30 labeling reactions of 20 µl each (35 % DIG-11-UTP substitution, 1 mM ATP, GTP, CTP, 0.65 mM UTP, 0.35 mM DIG-11-UTP).

#### Preparation of working solutions:

##### A. Preparation of 10 mM ATP/CTP/GTP working solution

- Thaw 100 mM ATP, 100 mM CTP and 100 mM GTP solutions on ice, vortex and spin-down briefly.
- Prepare a 1:10 dilution with Nuclease-Free Water to achieve a final concentration of 10 mM (e.g. 5 µl 100 mM ATP + 5 µl 100 mM CTP + 5 µl 100 mM GTP + 35 µl Nuclease-Free Water).
- 10 mM ATP/CTP/GTP working solution can be stored at -20°C. Avoid freeze/thaw cycles.

##### B. Preparation of 10 mM UTP working solution

- Thaw 100 mM UTP solution on ice, vortex and spin-down briefly.
- Prepare a 1:10 dilution with Nuclease-Free Water to achieve a final concentration of 10 mM (e.g. 5 µl 100 mM UTP + 45 µl Nuclease-Free Water).
- 10 mM UTP working solution can be stored at -20 °C. Avoid freeze/thaw cycles.

#### Transcription Reaction Setup (20µl):

1. Place T7 RNA Polymerase Mix on ice.
2. Thaw all remaining components at room temperature (RT), mix and pulse-spin in microfuge to collect solutions to bottom of tubes.
3. Assemble the reaction at room temperature in the following order:

Component	Volume	Final conc.
1. Nuclease-Free Water	X µl	



## DATA SHEET

Version: 02

Revision date: 13/03/2023

2.	T7 10X Reaction Buffer (10x)	2 $\mu$ l	1x
3.	DTT (100 mM)	2 $\mu$ l	10 mM
4.	10 mM ATP/CTP/GTP	2 $\mu$ l	1 mM
5.	10 mM UTP	1.3 $\mu$ l	0.65 mM
6.	10 mM DIG-11-UTP	0.7 $\mu$ l	0.35 mM
7.	Template DNA	X $\mu$ l	0.5-1 $\mu$ g
8.	T7 RNA Polymerase Mix	2 $\mu$ l	

4. Mix thoroughly, pulse-spin in microfuge. Incubate at 37°C for 30 minutes in the dark.

Depending on the RNA probe individual optimization may increase product yield (2h – 4h at 37 °C).

Note that, modified ribonucleotides reduce transcription efficiency; therefore, lower transcription yields should be expected as compared to transcription using unmodified NTP.

5. **Optional:** Once the in vitro transcription reaction is finished, treat the sample with RNase-Free DNase I to remove template DNA and purify the RNA. *For many applications it may not be necessary.*

Proceed with a DNase digestion as following:

Add RNase-free DNase I (not provided) plus the respective reaction buffer directly to the transcription reaction; mix and incubate for 15 minutes at 37 °C. Stop DNase activity by adding 1  $\mu$ L of EDTA at 0.5 M.

6. Proceed with purification of synthesized RNA or analysis of transcription products by gel electrophoresis.

