

DATA SHEET


Version: 2
Revision date: 24/02/2023Tlf: +34 983 54 85 63
info@canvaxbiotech.comwww.canvaxbiotech.com

1. Identification

Product name	Puromycin dihydrochloride Solution
	10ml
Cat. No	AB027

2. Description

Puromycin Hydrochloride is the hydrochloride salt form of puromycin, an aminonuclease antibiotic agent isolated from the bacterium *Streptomyces alboniger*. Toxic to both eukaryotic and prokaryotic cells. Acting as an analog of the 3' terminal end of aminoacyl-tRNA, puromycin incorporates itself into a growing polypeptide chain and causes its premature termination, thereby inhibiting protein synthesis. Prevents growth of bacteria, protozoa, algae, and mammalian cells.



3. Specifications

- CAS Number: 58-58-2
- Chemical Formula: C₂₂H₂₉N₇O₅·2HCl
- Molecular Weight: 544.44
- Purity (HPLC): > 98%

4. Storage specifications

Store at -20° C

5. Applications

- Puromycin is widely used in cell biology as a selection antibiotic agent in mammalian cell culture systems.
- Useful in CRISPR/Cas9 mammalian gene editing by selecting for successful Cas9-induced knock-in with puromycin resistance gene.
-

6. Further information

Product Use Limitation This product is developed, designed and sold exclusively only for research purposes use. The product was not tested for use in diagnostics or for drug development, nor is it suitable for administration to humans or animals.

Disclaimer The information provided in this Data Sheet is correct to the best of our knowledge and belief at the date of publication. This information is intended only as a guide and should not be taken as a warranty or quality specification. Canvax Reagents S.L.U. shall not be held liable for any damage resulting from handling or from contact with the above product.

