

DATA SHEET

Version: 2
Revision date: 17/02/2023

Canvax Reagents, S.L.U.
Luis de Mercado Street, 19
Boecillo Technological Park
47151, Boecillo
Valladolid, Spain.

Tlf: +34 983 54 85 63
info@canvaxbiotech.com

www.canvaxbiotech.com

1. Identification

Product name	Phosphate Buffered with Tween™ 20, pH 7.4 25L (50 tablets)
Cat. No	BR0007

2. Description

Phosphate-buffered Saline with Tween™-20 (PBS-T) is a water-based salt solution ideal for use in sample preparation and as a wash buffer in general immunoassay applications. It contains Tween 20, a non-ionic detergent additive that reduces non-specific binding and protein-protein interaction during the wash step in protein and immunoassay procedures such as ELISA and Western blotting. Decreasing the non-specific binding and staining makes ELISA results and blots easier to interpret.

3. Specifications

Chemicals: Analytical grade.

Format: Exactly pre-weighed tablets

Composition: 0.01M phosphate buffered saline, 0.14 M NaCl, 0.0027 M KCl, and 0.05% Tween 20
pH 7.4 ± 0,2 at 25°C.

4. Shipping and Storage specifications

PBS-T buffer is shipped at room temperature. Store the product in a dry place at room temperature.

5. Applications

- Wash buffer in immunolabelling techniques, such as ELISA and Western blotting
- ELISA and Western blotting
- Blocking buffer for plate-based assays
- Protein-plate coating

6. Directions for use

Dissolve one tablet in 500mL of water.

7. Further information

Product Use Limitation This product is developed, designed, and sold exclusively only for research purposes use. The product was not tested for use in diagnostics or for drug development, nor is it suitable for administration to humans or animals.

Disclaimer The information provided in this Data Sheet is correct to the best of our knowledge and belief at the date of publication. This information is intended only as a guide and should not be taken as a warranty or quality specification. Canvax Reagents S.L.U. shall not be held liable for any damage resulting from handling or from contact with the above product.

