

DATA SHEET

Version: 2
Revision date: 17/02/23

1. Identification

Product name	3 M Sodium Acetate (pH 5.2)
	100 mL
Cat. No	BR0080

2. Description

3M Sodium acetate solution, pH adjusted to 5.2 with glacial acetic acid, Sterile, 100mL.

3. Specifications

Chemicals: Biotechnology grade
Format: Solution
Concentration: 3M
Volume: 100 mL
pH: 5.2 ± 0.05 at 25 °C.

4. Shipping and Storage specifications

Sodium Acetate is shipped at room temperature. Store the product at Room Temperature.

5. Applications

Ideal for a wide range of molecular biology applications.

6. Further information

Product This product is developed, designed and sold exclusively only for research purposes use.
Use The product was not tested for use in diagnostics or for drug development, nor is it suitable
Limitations for administration to humans or animals.
Disclaimer The information provided in this Data Sheet is correct to the best of our knowledge and belief at the date of publication. This information is intended only as a guide and should not be taken as a warranty or quality specification. Canvax Reagents S.L.U. shall not be held liable for any damage resulting from handling or from contact with the above product.

