

DATA SHEET

Version: 1
Revision date: 02/06/2025

1. Identification

Product name 5M Sodium Chloride Solution, Sterile**Cat. No** BR0397 1 L
BR0398 5 L

2. Description

Enhance your processes with **5M Sodium Chloride Solution**, a high-purity, sterile sodium chloride solution designed to ensure superior quality and performance in scientific and medical applications. Manufactured under stringent quality controls, this solution provides stability and safety in research, laboratory, and production environments.

3. Specifications

Format: Liquid, Ready-to-use**Composition:** DI Water (70.78%) and Sodium Chloride (29.22%)**Sterilization:** purified using a 0.22 µm filtration process [aseptically filled]

4. Shipping and Storage specifications

This product is shipped at Ambient Temperature. Store the product in a dry place at RT (15 – 25° C).

5. Applications

Designed for the most commonly used **Research and Biomanufacturing applications**, such as:

- Cell Culture & Microbiology
- Buffer Preparation
- Biomanufacturing
- Industrial applications

6. Quality

Manufactured in an **ISO 13485:2016 and GMP-certified facility**, located in Boecillo, Valladolid (Spain), under ISO9001:2016 procedures to ensure the highest quality standards.

7. Further information

Product Use Limitation This product is developed, designed and sold exclusively only for research purposes use. The product was not tested for use in diagnostics or for drug development, nor is it suitable for administration to humans or animals.

Disclaimer The information provided in this Data Sheet is correct to the best of our knowledge and belief at the date of publication. This information is intended only as a guide and should not be taken as a warranty or quality specification. Canvax Reagents S.L.U. shall not be held liable for any damage resulting from handling or from contact with the above product.

