

Data Sheet

Version: v1
Revision date: 06/05/2024

1. Identification

| | |
|---------------------------|----------------------------|
| Product name | WFI Quality Water |
| Manufacturing Site | Boecillo; Valladolid Spain |

2. Description

Water for Injection (WFI) Quality Water, Animal-derived component-free, meets the most stringent USP and/or EP compendial grade specifications. Manufactured in an **ISO 13485:2016 and GMP-certified facility**, located in Boecillo, Valladolid (Spain), following GMP compliant procedures to ensure the highest quality standards. It provides you the confidence to achieve exceptional results, with unparalleled purity and reliability with a secure supply chain.

3. Size, SKU and Packaging

| | |
|--------------|---|
| E0335 | 1,000 mL served in a PET Bottle. |
| E0361 | 20 L served in a graduated square container. |
| E0371 | 200 L served in a single-use bag with plastic drum. |

4. Storage & Stability specifications

- **Shipment** at Ambient Temperature.
- **Storage:** Room Temperature (15 - 30 °C).
- **Shelf life:** 24 months from data of manufacturing, based on our Accelerated Stability Testing Program, with both Normal and Accelerated Stability Testing Programs ongoing.

5. Applications

For Research Use Only or Further Manufacturing. Not for diagnostic use or direct use in humans or animals.

6. Quality

Manufactured in an **ISO 13485:2016 and GMP-certified facility**, located in Boecillo, Valladolid (Spain), under ISO9001:2016 procedures to ensure the highest quality standards. Our processes include:

- Batch-to-batch consistency.
- Facilities audited regularly Quality audits.
- Lot specific CoAs and QA review previously to release a product.
- Material traceability, documented procedures and documented employee training.
- Equipment maintenance and monitoring records.



Data Sheet

Version: v1
Revision date: 06/05/2024

Canvax Reagents, S.L.U.
Luis de Mercado Street, 19
Boecillo Technological Park
47151, Boecillo Valladolid, Spain.
Tlf: +34 983 54 85 63
info@canvaxbiotech.com
www.canvaxbiotech.com

7. Further information

- Product Use Limitation** This product is developed, designed and sold exclusively only for research purposes use and Further Manufacturing. The product was not tested for use in diagnostics or for drug development, nor is it suitable for administration to humans or animals.
- Disclaimer** The information provided in this Data Sheet is correct to the best of our knowledge and belief at the date of publication. This information is intended only as a guide and should not be taken as a warranty or quality specification. Canvax Reagents S.L.U. shall not be held liable for any damage resulting from handling or from contact with the above product. It is the end user's responsibility to ensure that the final product meets the requirements of the application for which it is to be used.

