

DATA SHEET

Version: 2
Revision date: 10/04/2023

1. Identification

Product name	6X B/Loading Buffer
	5 x 1 mL
Cat. No	L0031

2. Description

6X B/Loading Buffer is used as a loading dye for visual tracking of DNA migration during electrophoresis. It contains Bromophenol blue. Bromophenol blue migrates fast in the agarose gel and corresponds to the migration of 300 - 500 bp long DNA fragment in 1% agarose gel. The EDTA is included in the solution to protect sample from nuclease degradation. Glycerol is added to provide high density to solution. Due to its high density, sample settled at the bottom of the well. It also helps DNA sample to be confined in the well without diffusing out from the well.



3. Concentration

0.03% bromophenol blue, 60% Glycerol, 10 mM Tris-HCl (pH 7.6) and 60 mM EDTA.

4. Storage

Store the solution at -20°C for long time. Solution can be stored at 2 - 8 °C for few weeks.

5. Further information

Product Use Limitations	This product is developed, designed and sold exclusively only for research purposes use. The product was not tested for use in diagnostics or for drug development, nor is it suitable for administration to humans or animals.
Disclaimer	The information provided in this Data Sheet is correct to the best of our knowledge and belief at the date of publication. This information is intended only as a guide and should not be taken as a warranty or quality specification. Canvax Reagents S.L.U. shall not be held liable for any damage resulting from handling or from contact with the above product.

