

DATA SHEET

Version: 2
Revision date: 10/03/23www.canvaxbiotech.com

1. Identification

Product name **TruePure™ dATP**
1 mL, 100 mM

Cat. No N0093-XL

2. Description

Ultrapure dATP supplied as clear aqueous solution (pH 8.5).

3. Specifications

Formula: C₁₀H₁₆N₅O₁₂P₃
Molecular Weight: 491.18 g/mol
λ_{max} pH 7.0 = 259 nm
ε at λ_{max}, pH 7.0 = 15.4 mmol⁻¹ cm⁻¹
Purity: > 99% confirmed by HPLC
Quality Control Specifications:

- 18 kb long range PCR (template dilution series): suitable
- Contamination with bacterial and human DNA: not detectable
- Activity of DNase, Protease or Phosphatase: not detectable

4. Storage specifications

Store at -20°C, short term (up to one week) exposure to ambient temperature is possible.

5. Applications

For use in all molecular biology applications, including PCR, real-time PCR, high fidelity and long PCR, cDNA synthesis, RT-PCR, and DNA sequencing.

The synthesis of a complementary DNA/cDNA strand by enzymatic polymerization of deoxynucleotide triphosphates (dNTPs) is the basic principle of DNA amplification, sequencing and cDNA synthesis techniques.

6. Further information

Product Use Limitations This product is developed, designed and sold exclusively only for research purposes use. The product was not tested for use in diagnostics or for drug development, nor is it suitable for administration to humans or animals.

Disclaimer The information provided in this Data Sheet is correct to the best of our knowledge and belief at the date of publication. This information is intended only as a guide and should not be taken as a warranty or quality specification. Canvax Reagents S.L.U. shall not be held liable for any damage resulting from handling or from contact with the above product.

